



SPECIFICATION

Product Code: 31390D



Product Name: BUFFER CAPSULES pH 4.0

Direction for use:

Open the capsule by twisting and pulling apart.

Dissolve the contents in 10 ml of CO₂ – free water
and then make it with CO₂ – free water upto 100 ml.

This solution will have

pH of solution (20⁰C) 4.0 ± 0.05

MOLYCHEM Office: 78/80 Babu Genu Road, Souri Bldg., 2nd Floor, Mumbai – 400 002 (INDIA)

Tel: +91 22 4066 8111, 2205 7772 . Fax: + 91 22 2203 4442. Email : info@molychem.net

Works: Plot No.B/9, Badlapur M.I.D.C.,OPP.Shobhavat Textile Co.,

Kulgaon, Badlapur (E), Dist.-Thane-421 503